

# INFECTION CONTROL Why Buy ?

When you can Hire the Instruments and the expertise

## Endoscope Prewash



**Let the Wash Wizard take the strain – not your hands !!**

- ❑ Carries out Flush /Rinse Pre-Wash procedures from pre-prepared detergents or concentrates.
- ❑ Follows Olympus Guidelines \* on Pre-wash flush and rinse procedures
- ❑ Does not use Syringes for the main lumens
- ❑ Carries out the whole air water lumen FLUSH process within 6 minutes

- ❑ Will not over-pressurise the scopes [unlike the syringe flush process!!]
- ❑ Simple Mimic screen tells of flush programme status
- ❑ SUPROX Technology for Self disinfect as standard

### Ask about our free training workshops

\* Olympus Guideline- Sequence

90cc detergent / 90cc water washes followed by an air flush – air and water channels  
 30cc detergent wash / 30cc water washes followed by an air flush – Elevator Channel

- ❑ **Options**
- ❑ Self contained air pump
- ❑ Eductor . 'Dosatron' for auto dilution of Detergent
- ❑ Silver ion residual sanitising

### Rent or Buy

Ask For Details

**PASA**

[www.medipureireland.net](http://www.medipureireland.net)

Medipure Ireland Ltd Head Office  
 Hickeys Yard , Old Ashbourne Rd The ward Dublin 11  
 Web [www.suprox.com](http://www.suprox.com) [www.medipure.co.uk](http://www.medipure.co.uk)  
 Company Number - 3996456 FS80797  
 Tel 0870 203 2030 01244 332441 / 332887

**CE**

0086

